

Insight on Color

Vol. 18, No. 5

Accessing EasyMatch QC Data in External Spreadsheets and Databases

Data obtained using a HunterLab sensor and EasyMatch QC software can be used in many ways. In addition to saving data, manipulating displays, and printing while within EasyMatch QC, you can also export or transfer the data to other programs for manipulation. Several methods for doing this are described below.

Copying and Pasting to Windows Spreadsheets

Data contained in EasyMatch QC's Color Data Table and Spectral Data Table views are already in spreadsheet format. You may copy the cells desired in EasyMatch QC and then paste them into other spreadsheet programs.

- 1. Configure the Color Data Table or Spectral Data Table to show the parameters, standards, and samples desired.
- 2. Highlight the EasyMatch QC data for copying by clicking and holding the left mouse button down while dragging the cursor across the cells desired.





3. To put the cells on the Windows clipboard, click the right mouse button in the spreadsheet view and select **Copy to Clipboard** from the menu that appears.

EasyMatchQC - [Untitled File Edit View Measuremer	Job2] [DataBase:EZQC Its Options Sensor Window	2				
	🚽 🛃 🛃 🚱		2			
+ Standard Empty		L×	a×	b*	dE×	_
	+Tolerances	0.00	0.00	0.00	0.00	
- Standard 1	-Tolerances	0.00	0.00	0.00	0.00	
- 🤹 Samples	Sample 2	81.87	-32.81	59.13	0.08	
	Sample 3	81.86	-32.77	59.29	0.16	
Sample 2	Sample 4	81.90	-32,93	59.13	0.20	
🧾 Sample 3	Sample 5	81.90	-32	Copy to Clipboard		
	Sample 2	81.87	-32	Configure	-10	
- 🈏 Sample 4	Sample 3	81.86	-32	Coningare		
- 🌅 Sample 5	Sample 4	81.90	-32	Print Preview		
20/20/20	Sample 5	81.90	-32	Print this view		-
Sample 2	H + + H D65/10 / F02	2/10 / A/10 /		12 32		
🧭 Sample 3	Color Data Table			Change View	•	
The second se				Delete View		
Sample 4	100 % R/T 90		of the local division of the local divisiono	Split View Vertically ar	1.1.1.1.1.1.1.1.1.1.1.1.1.1.1.1.1.1.1.	
🧭 Sample 5	Standard 1 80		/	Split View Horizontally	and Add	
	Mn/Mx : 11.52(450) / { 60 Sample 2 40		/	Maximize this View		
	Mn/Mx : 11.55(450) / E 20 Sample 3 10 400		500	Navelenoth (nm)	00	700
	Spectral Plot(Reflectance/Trans	mittance)				
🔶 Untitled Job2						
Ready Current Sens	or :LabScan XE "LSXEDEMO"	Current Stdz.N	lode : Mode	#1 - Reflectance - 3.	175 mm - Nominal	- Stdz. Expires 🏑

4. Open your Windows spreadsheet, click on any cell, and select **Paste** from the **Edit** menu to paste the data into the spreadsheet. The example given here uses Microsoft Excel 2003.

🔊 🖌 Eile	Edit	View Insert Form		la Data	110 malaura	Help A			Turno o gi	estion for help		
	-			_	_		-			-		r ^
1 🞽		Can't Undo Ctrl+	z X	la 🖪	- 🏈 🛛 🖻	7 - (2 - 1	🔒 Σ ᠇		100%	• 🕐 🚽		
🔁 🕲 '	U	Can't <u>R</u> epeat Ctrl+	Y II	12 180	Reply with	Changes	End Review.					
Arial	X	Cu <u>t</u> Ctrl+	× n		-	\$ %	€.0 .00		🛛 • 🖄 • 🛓	A		
A1		Copy Ctrl+	c 🎽			Ψ /0 /	.00 -+.0		<u> </u>	- 5		
	-	Office Clip <u>b</u> oard	[E	F	G	Н		J	К	-
1	2	Paste Ctrl+			-						1X	_^
2		Paste Special										
3		-										
4		Paste as <u>H</u> yperlink	_									
6		Fill	+									
7		Cle <u>a</u> r	•									
8		Delete										
9		Delete Sheet										
10 11		Move or Copy Sheet	-									_
10	æ	Eind Ctrl+	-									
13	ara	-										
14		Replace Ctrl+I										
15		Go To Ctrl+	5	-								
16 17		Lin <u>k</u> s	-									_
18		Object		_								
19	-											~
	\ sł	neet1 / Sheet2 / She	et3 /				<				>	

5. Manipulate the data and/or perform calculations as desired.



	<u>Eile E</u> dit	<u>V</u> iew <u>I</u> nse	rt F <u>o</u> rmat	Tools D	ata <u>W</u> indov	v <u>H</u> elp A	do <u>b</u> e PDF		Type a c	uestion for help	d	8 ×
	📬 🖬 🗋	a a a i	ية 💝 ا 🔊	1 X D	2 🖋 🛛	•) • (•! •	😣 Σ 🔹		100%	• 💿 📮		
	* * 2	1 🗞 🏹 🛙	5 81 🕽		♥ Reply with	h Changes	End Review.					
Aria	l	- 10	- B I	U 🔳	= = •	\$ %	, .0 .00 0.€ 00.		🖽 • 🖄 • ,	A -		
·	A7		s.									
	A	В	С	D	E	F	G	Н	1	J	K	-
1	81.87	-32.81	59.13	0.08								
2	81.86	-32.77	59.29	0.16								
3	81.90	-32.93	59.13	0.20								
4	81.90	-32.76	59.14	0.05								
5	81.87	-32.81	59.13	0.08								
6	81.88	-32.82	59.16	0.11	Average							
7												
8		8	1									
9												
10												
11												
12												
13												
14		1										
15												t l
16	- C	1										t l
17												
18												
19												-
	▶ N\Sh	eet1 / Shee	t2 / Sheet	3/			<				>	
Ready			,,				1000					-7.4

Opening EasyMatch QC's Database Files Using Another Database or Spreadsheet Program

EasyMatch QC gives you the option of saving to either a Microsoft Access or SQL Server database, and that database may be opened directly through Access/SQL in addition to through EasyMatch QC. The table titled "tSensorTypeSampleSpecData" contains the raw spectral data for each sample stored. "tSensorTypeStandardSpecData" contains the same information for each standard stored. It is this information that is then converted to color scale and index information as needed for display on your screen using formulas internal to EasyMatch QC. Only this spectral data is directly available in the database, so if you need to access color scale or index information, you will need to either copy and paste as described above or export this information as described below.

Note: For security reasons, the database may not be directly opened in Access or SQL if you are running EasyMatch QC-ER.

ASCII Export

You may directly export data from EasyMatch QC to an ASCII file that may then be opened in other programs, such as spreadsheet and database programs. The **ASCII Export/Configure** command in the **Options** menu allows you to select any color scale, color difference scale, index, spectral data, etc. information desired for export to an ASCII file. **ASCII Export/Export Data** in the **Options** menu is then used to export the chosen parameters for the standards and samples desired to an ASCII text file. Automatic export of measurements as they are made is also available. Complete instructions for this procedure are included in your EasyMatch QC User's Manual.



elected Items			Scale	Illuminant/Observer
TEXT\Measurement ID		<<	CIELAB	D65/10 💌
SCALE\CIELAB\D65/10 SPECTRAL\10 nm	*		Delta Scale	MI Illuminant
		<<	dL*	A/2 💌
	×		Indices	
			457nm Brightnes: 💌	
		>>>	Spectral Data	
		<<	10 nm 💌	
			Text Fields	
		<<	Measurement ID	-
			Procedure	
		<<	Opacity 💌	
			Delimiter	
			C Comma	
			Tab Other	ASCII Code
Auto Export New Data		ок	144	
Auto Export Uses This File:			Quotes Around 1	ext Fields
C:\EasyMatchQCAsciiExport.txt		Cancel	🔽 Include Column I	

Note: If you select a Tab Delimiter and use an extension such as .XLS or .CSV, your exported file may be opened directly into Microsoft Excel. If you select a Comma Delimiter and use an extension such as .CSV, your exported file may be imported into a new table in Microsoft Access 2000 or higher.

Data Send

You may send data from EasyMatch QC directly to a communications port or IP address so long as software is available at that address to accept the data. The **Data Send/Configure** command in the **Options** menu allows you to select any color scale, color difference scale, index, spectral data, etc. information desired for sending. **Data Send/Send Data** in the **Options** menu is then used to send the chosen parameters for the standards and samples desired to the configured address. Automatic sending of measurements as they are made is also available. Complete instructions for this procedure are included in your EasyMatch QC User's Manual.





For Additional Information Contact:

Technical Services Department Hunter Associates Laboratory, Inc. 11491 Sunset Hills Road Reston, Virginia 20190 Telephone: 703-471-6870 FAX: 703-471-4237 www.hunterlab.com

