Method #1: For Direct Connection between Aeros and Computer with Ethernet Cable (No Network required)

Method #2: For Direct Connection between Aeros and Computer with RS-232 Ethernet Adapter (No Network required)

Method #3: Through Network (Ethernet cable or WiFi)

# Method #1: For Direct Connection between Aeros and Computer with Ethernet Cable (No network required)

- 1) Materials Needed
  - Hardware needed: Ethernet cable, and USB Ethernet adapter if you want to have multiple Ethernet cables connecting to the computer.



Figure 1- Ethernet Cable



Figure 2- USB Ethernet Adapter

2) Connect Aeros to Computer:

 Plug Ethernet cable (Figure 1) into RJ-45 Ethernet connection at rear of Aeros. Plug other end of cable into the computer directly or through USB Ethernet Adapter (Figure 2).



Figure 3- Rear View of Aeros

- 3) Configure IP address of the Aeros: (Note: Requires Essentials Rev 14 or higher)
  - Check IP address in PC: Open Command Prompt in the PC. Type in "ipconfig" and check the "autoconfiguration IPv4 Address" as well as the "Subnet Mask".




Figure 4- Computer IP address

 Configure IP address in Aeros: Open Aeros Essentials, go to Preferences/Configure Network Settings. Uncheck "Use DHCP for Ethernet Config". Type in IP address and Subnet Mask manually, then apply. The IP address here should be exact same as the "autoconfiguration IPv4 Address" in the PC, except changing the last two digits. Restart Aeros to get network setting applied.





Figure 5- Assign IP address to Aeros

- Press Apply on the Ethernet Configuration
  - Check "Automatic Network Data Export" and click to configure. For a direct connection between Aeros and computer with Ethernet cable, set up the Aeros as a client. Port Number is 10001 (Or 11111, Just make sure that you use exact same port number in Essentials and PC). Enter the PC IP



# 4) Configure the Computer Using HyperTerminal Software

- Go to the computer and open HyperTerminal.
- Make a new connection:
  - 1. Enter a name for the connection
  - 2. Connect using TCP/IP(Winsock)
  - 3. Enter the PC IP address and Port Number 10001 (Or 11111, Just make sure that you use exact same port number in Essentials and PC)

#### 5) Send Data from the Aeros:

- Configure the Aeros for the Color Data Screen:
  - 1. Select Color Scales, Indices & Illuminant/Obs (Workspace>Color Scales)
  - 2. Standardize the instrument. (Workspace>Standardize)
- Read sample and view the data on the computer.




Figure 7- Data from Aeros through Ethernet Connection

The data string is shown as follows:

<STX><PACKET SIZE>\$T<SCALE-LABEL>\$T<iII/Obs>\$T<LABEL
NAME1><\$T><VALUE1>\$T<LABEL NAME2><\$T><VALUE2><\$T><LABEL
NAME3>\$T<VALUE3>\$T ..... <LABEL NAME N><\$T>

<VALUE N>\$T<CHKSUM><ETX>

Where,

<STX> is the Start of Text (value =0x02)

 $\langle ETX \rangle$  is the End of Text (value =0x03)

\$T is the default delimiter.

<SCALE-LABEL> is the Scale Label (e.g. CIELAB)

<ill/Obs> is the ill/Obs name (e.g. D65/10)

<PACKET SIZE> is the Total size (HEX String) of the Packet excluding the <STX> and <ETX>

<LABEL NAME> is the label name (e.g. L\*, a\*, b\*, dE\* etc..)

<VALUE> is the value for the preceding Label Name

<CHKSUM> is the Checksum (HEX String) - the sum of all the ASCII values in the total packet play load starting from <PACKET-SIZE> and till <CHKSUM>

# Method #2 - For Direct Connection between Aeros and Computer with RS-232 Ethernet Adapter

# 1) Materials Needed

 Hardware needed: Ethernet cable, Crossover adapter, Ethernet to RS-232 Connector, RS-232 to USB (optional)



Figure 8- Ethernet Cable



Figure 9- Crossover Adapter



Figure 10- Ethernet to RS-232 converter for Connection via Serial port



Figure 11- RS-232 to USB converter for connection via USB port

# 2) Configure Ethernet to RS-232

- Set up Ethernet to the RS-232 adaptor with static IP address and Port Number.
- In this example, the adapter was configured with IP address of 192.168.0.100 and port 10001 (Or 11111, Just make sure that you use exact same port number in Essentials and PC).

#### 3) Connect Aeros to Computer:

• Plug Ethernet cable (Figure 8) into RJ-45 Ethernet connection at rear of Aeros. Plug other end of cable into Crossover Adapter (Figure 9).



Figure 12- Rear View of Aeros

- Plug Crossover Adapter into Ethernet port of Ethernet to RS-232 adapter (Figure 10).
- Plug Ethernet to RS-232 adapter into serial port of computer (if D-9 serial port is present) or into RS-232 to USB converter for connection to USB port. Plug power into Ethernet to RS-232 adapter.




Figure 13- Cable Configuration for Direct Computer Connection

- 4) Configure the Aeros: (Requires Essentials Rev 14 or higher)
  - Configure the Ethernet port of Aeros. Select Jobs Menu> Preferences and Select 'Configure Network Settings'. Uncheck "Use DHCP for Ethernet Config" and enter a valid IP address for the Ethernet port. In this example, the following parameters are selected.




Figure 14- Configuration Parameters for Ethernet

- Press Apply on the Ethernet Configuration and then Apply on the Preferences Page to complete.
- Go to Jobs Menu> Preferences, check "Automatic Network Data Export" and click to configure.

Figure 15- Preferences>Automatic Network Data Export

- For a direct connection between Aeros and data collection computer, set up the Aeros as a Client.
- Set the IP Address to match the settings of the Ethernet to RS-232 converter, in this case 192.168.0.100 and the Port as 10001 (Or 11111, Just make sure that you use exact same port number in Essentials and PC).
- Press Apply and then press Apply on the Read Options screen to continue.

Figure 16- Read Options Export

## 5) Configure the Computer:

- Connection configurations differ depending on data collection software. In this example, Hyperterminal is used to demonstrate connectivity.
- The data collection computer will be set up as a Server.
- Connect as follows:
  - 1. Select the Comport that represents the USB or Serial port connection.
  - 2. Define the Connection:
  - 3. Bits per second: 9600
  - 4. Data Bits: 8
  - 5. Parity: None
  - 6. Stop Bits: 1
  - 7. Flow Control: None

## 6) Send Data from the Aeros:

- Configure the Color Data Table with the color scale and parameters to be transmitted.
  - 1. Standardize the Instrument
  - 2. Select Color Scales, Indices & Illuminant/Obs (Workspace>Color Scales)
  - 3. Hit Read Button, data will be transferred to the computer.




Figure 17- Data Output

The data string is shown as follows:

<STX><PACKET SIZE>\$T<SCALE-LABEL>\$T<iII/Obs>\$T<LABEL
NAME1><\$T><VALUE1>\$T<LABEL NAME2><\$T><VALUE2><\$T><LABEL
NAME3>\$T

<VALUE N>\$T<CHKSUM><ETX>

Where.

<STX> is the Start of Text (value =0x02)

<ETX> is the End of Text (value =0x03)

\$T is the default delimiter.

<SCALE-LABEL> is the Scale Label (e.g. CIELAB)

- <ill/Obs> is the ill/Obs name (e.g. D65/10)
- <PACKET SIZE> is the Total size (HEX String) of the Packet excluding the <STX> and <ETX>
- <LABEL NAME> is the label name (e.g. L\*, a\*, b\*, dE\* etc..)
- <VALUE> is the value for the preceding Label Name
- <CHKSUM> is the Checksum (HEX String) the sum of all the ASCII values in the total packet play load starting from <PACKET-SIZE> and till <CHKSUM>

# Method #3 – To export data through network

Connect Aeros to a Network. You can connect Aeros to a network hub using the Ethernet cable or connect Aeros to a network hotspot through a WiFi connection. The computer must be connected to the same network as the Aeros.

## Option A: Connect Aeros to the network with Ethernet cable.

- To connect Aeros to network, go to Workspaces> Preferences and Select 'Config Network Settings'. Select Ethernet configuration and check "Use DHCP for Ethernet Config". Please write down the IP address showing in the Ethernet Setting dialog. You can also check the IP address of Aeros in Jobs/About/Info.
- 2. Go back to Workspaces/Preferences and select *Auto Network Data Export* Measurement using a check and select "Config" button. Choose Aeros as *Server* and *Port number as 11111* (Or 10001, Just make sure that you use exact same port number in Essentials and PC).
- 3. You can also choose a delimiter to mark your data.
- 4. Configure the terminal program such as Hyperterminal in the computer with the following settings:
- Set computer as client
- Enter the IP address of Aeros which is recorded in step 1 above
- Put the port number as "11111" (Or 10001, Just make sure that you use exact same port number in Essentials and PC).

After all the setting, you should be ready to get the data exported from Aeros to the computer.

# Optional B: Connect to a hotspot through the WiFi connection

Note: This option only works for ARS00016 and above. All Aeros below this serial number does not have WiFi ability built-in.

- To connect Aeros to network, go to Workspaces> Preferences and select Config Network Settings. Select Configure WiFi Settings and the WiFi configuration dialog will be prompted. Please search and connect to the available WiFi and write down the IP address showing in this dialog. After the WiFi configuration, please click the floating Back Button to go back to Essentials app.
- Go back to Workspaces/Preferences and select Auto Network Data
   Export Measurement using a check and select "Config" button. Choose Aeros
   as Server and Port number as 11111 (Or 10001, Just make sure that you use exact
   same port number in Essentials and PC). You can also choose a delimiter to mark your
   data.
- 3. Configure the terminal program such as Hyperterminal in the computer with the following settings:
- Set computer as client

- Enter the IP address of Aeros which is recorded in step 1 above
  Put the port number as "11111" (Or 10001, Just make sure that you use exact same port number in Essentials and PC).

After all the setting, you should be ready to get the data exported from Aeros to the computer.