

# CFEZ/MSEZ Instrument Connecting and Disconnecting

---

**FAQ: "Help! I recently ran updates on my computer and now my instrument keeps connecting and disconnecting."**

This usually has something to do with the Enhanced Power Management setting in Windows 8. Please read the blog note below for instructions on how to disable this feature.

<http://blogs.msdn.com/b/usbcareblog/archive/2013/11/08/help-after-installing-windows-8-1-my-usb-device-doesn-t-charge-or-it-disconnects-and-reconnects-frequently.aspx>

A power saving workaround was implemented in EasyMatch QC 4.78 and higher to deal with this automatically.