If the "Login Required" box is checked under Options/Application Preferences/Startup, the user is prompted to log into EasyMatch QC.

As noted in the EasyMatch QC User's Reference Manual, the default EZMQC login and password with full admin privileges to operate all features within EZMQC are:

User Name: admin (not case sensitive)

Password: **hunterlab** (all small letters, case sensitive)

This default login and password can be removed or changed by the user.

If this has been done, the **backdoor login and password** can be used for access and to reset login/password under Options/System Configuration/User Manager.

User Name: hunterlab (not case sensitive)

Password: **marigold** (all small letters, case sensitive)

This backdoor login and password has full admin privileges to use all features within EasyMatch QC and cannot be removed or changed by the user. You can create new admin users in EasyMatch QC when using this backdoor entry.